
ELI 380 ELECTROCARDIOGRAPH EL CTROCARDIOGRAPH ELI 380 ELECTRO CARDIOGRAPH ELI 380 ELECTROCA RDIOGRAPH ELI 380 ELECTROCA OGRAPH ELI 380 ELECTROCARDIC



ELI 380 ELECTROCARDIOGRAPH

ELI 380 Solving the challenges of TODAY.



KEY ADVANTAGES:

CLINICAL SUPERIORITY


- The world-renowned VERITAS[™] Resting ECG Interpretation Algorithm offers impressive diagnostic ECG performance.
- Best-in-class signal acquisition technology exceeds published requirements providing increased ECG data accuracy.
- Mortara Critical Values enable a clinician to quickly and easily identify high-risk patients needing immediate attention.

STREAMLINED WORKFLOW

- One-touch "SYNC" function transmits ECGs and receives orders for specific care areas.
- Demographic pre-population via DICOM[®] Modality Worklist, ADT or barcode scanner reduces the steps necessary to get started.
- STAT ECG function enables ECG acquisition with the touch of a button.

SOLUTION-BASED DESIGN

- Large, bright 17-inch color display with backlight LED maximizes available data and visibility.
- Glass keyboard with touchpad promotes improved infection control compared to traditional keyboards.
- Lead reversal detection identifies limb and chest lead reversals in real time, a common error in high-volume, fast-paced departments.



ECG - AS EASY AS 1, 2, 3.

Alt A S D F G H J K L

1. HOOKUP

•

Using the WAM[™] (Wireless Acquisition Module), technicians no longer need to navigate the electrocardiograph within tight spaces. Instead, keep the cart outside of the care area and initiate the ECG right at the WAM, an ideal feature for isolation patients.

2. DEMOGRAPHIC IMPORT

• • •

Barcode-based workflow does away with time-consuming searches and makes upload of patient specific information fast and accurate.

-

3. ANALYSIS

Mortara's innovative Best 10 feature determines the most noise-free ECG data recorded and automatically displays the ECG, reducing both review time and the likelihood of repeat ECGs.

0



MORTARA INSTRUMENT, INC. 7865 North 86th Street

7865 North 86th Street Milwaukee, WI 53224 U.S.A. Tel: 414.354.1600 Tel: 800.231.7437 Service: 888.MORTARA Fax: 414.354.4760

Mortara

ELI™, Mortara®, VERITAS™, and WAM™ are trademarks or registered trademarks of Mortara Instrument, Inc.

DICOM[®] is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.

Ref 9513-189-00-ENG Rev 1

www.mortara.com

